

Safety data sheet
according to 1907/2006/EC, Article 31

Printing 27.01.2020

version number 5

Revision: 27.01.2020

SECTION 1: Identification of the substance/mixture and of the company/undertaking

***1.1 Product identifier**

***Trade name:** S3 KINEMATIC VISCOSITY,

***Article number:** REVISC-S3

***Registration number**

A registration number is not available for this substance as the substance or its uses are exempted from registration, the annual tonnage does not require a registration or the registration is envisaged for a later registration deadline.

***1.2 Relevant identified uses of the substance or mixture and uses advised against**

No further relevant information available.

***Application of the substance / the mixture** Laboratory Chemicals

***1.3 Details of the supplier of the safety data sheet**

***Manufacturer/Supplier:**

Reagecon Diagnostics Ltd.

Shannon Free Zone,

Shannon,

Co. Clare,

Ireland.

Tel +353 61 472622

Fax +353 61 472642

***Further information obtainable by contacting:** sds@reagecon.ie

***1.4 Emergency telephone number:**

For Hazardous Materials [or Dangerous Goods] Incident

Spill, Leak, Fire, Exposure, or Accident

Call CHEMTREC

For Ireland call +(353)-19014670

For Outside Ireland call +1 703-741-5970 / 1-800-424-9300 CCN849800

SECTION 2: Hazards identification

***2.1 Classification of the substance or mixture**

***Classification according to Regulation (EC) No 1272/2008**



GHS08 health hazard

Asp. Tox. 1 H304 May be fatal if swallowed and enters airways.

***2.2 Label elements**

***Labelling according to Regulation (EC) No 1272/2008**

The product is classified and labelled according to the CLP regulation.

***Hazard pictograms**



GHS08

***Signal word** Danger

***Hazard-determining components of labelling:**

Distillates (petroleum), hydrotreated light

dodecane

***Hazard statements**

H304 May be fatal if swallowed and enters airways.

***Precautionary statements**

P301+P310 IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

P331 Do NOT induce vomiting.

(Contd. on page 2)

Safety data sheet
according to 1907/2006/EC, Article 31

Printing 27.01.2020

version number 5

Revision: 27.01.2020



Trade name: S3 KINEMATIC VISCOSITY,

(Contd. of page 1)

P405 Store locked up.

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

2.3 Other hazards**Results of PBT and vPvB assessment*****PBT:** Not applicable.***vPvB:** Not applicable.**SECTION 3: Composition/information on ingredients*****3.2 Chemical characterisation: Mixtures*****Description:** Mixture of substances listed below with nonhazardous additions.***Dangerous components:**

EC number: 927-632-8	Distillates (petroleum), hydrotreated light	 Asp. Tox. 1, H304	50-100%
CAS: 112-40-3 EINECS: 203-967-9	dodecane	 Asp. Tox. 1, H304	10-25%

Additional information:** For the wording of the listed hazard phrases refer to section 16.**SECTION 4: First aid measures**4.1 Description of first aid measures*****After inhalation:**

Provide fresh air, warmth and rest. When breathing is difficult, properly trained personnel may assist affected person by administering oxygen. Obtain medical attention if any discomfort continues.

After skin contact:** Generally the product does not irritate the skin.After eye contact:**

Promptly wash eyes with plenty of water for up to 15 minutes. Open eyes wide apart and rinse well to remove any contact lenses. Do not remove contact lenses by hand. Continue to rinse. Get medical attention if symptoms persist.

***After swallowing:**

Do not induce vomiting; call for medical help immediately. Rinse mouth thoroughly with water and give large amounts of water to drink. Never give anything by mouth to an unconscious person.

Caution if victim vomits. Risk of aspiration. Keep airways free. Call a physician immediately.

4.2 Most important symptoms and effects, both acute and delayed** No further relevant information available.4.3 Indication of any immediate medical attention and special treatment needed**

No further relevant information available.

SECTION 5: Firefighting measures***5.1 Extinguishing media*****Suitable extinguishing agents:**

Indications shall be given whether any extinguishing media are inappropriate for a particular situation involving the substance or mixture

Use fire extinguishing methods suitable to surrounding conditions.

5.2 Special hazards arising from the substance or mixture** No further relevant information available.5.3 Advice for firefighters*****Protective equipment:**

Fire fighters should wear complete protective clothing including self-contained breathing apparatus. Do not breathe fumes. Keep containers cool by spraying with water if exposed to fire. Avoid run off to waterways and sewers

(Contd. on page 3)

Safety data sheet
according to 1907/2006/EC, Article 31

Printing 27.01.2020

version number 5

Revision: 27.01.2020

Trade name: S3 KINEMATIC VISCOSITY,

(Contd. of page 2)

SECTION 6: Accidental release measures***6.1 Personal precautions, protective equipment and emergency procedures**

Wear protective equipment as described in Section 8 below. Keep unprotected persons away.

6.2 Environmental precautions: Do not allow to enter sewers/ surface or ground water.**6.3 Methods and material for containment and cleaning up:**

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

Dispose contaminated material as waste according to item 13.

Ensure adequate ventilation.

***6.4 Reference to other sections**

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

SECTION 7: Handling and storage***7.1 Precautions for safe handling**

Ensure good ventilation/exhaustion at the workplace.

Prevent formation of aerosols.

Information about fire - and explosion protection: No special measures required.**7.2 Conditions for safe storage, including any incompatibilities*****Storage:*****Requirements to be met by storerooms and receptacles: No special requirements.*****Information about storage in one common storage facility: Not required.*****Further information about storage conditions: None.*****7.3 Specific end use(s) No further relevant information available.****SECTION 8: Exposure controls/personal protection*****8.1 Control parameters*****Ingredients with limit values that require monitoring at the workplace:**

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

Additional information: The lists valid during the making were used as basis.**8.2 Exposure controls*****Personal protective equipment:*****General protective and hygienic measures:**

Keep away from foodstuffs, beverages and feed.

Wash hands before breaks and at the end of work.

***Respiratory protection:**

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.

Where risk assessment shows air-purifying respirators are appropriate use a respirator with multi-purpose combination (US) or type ABEK (EN14387) respirator cartridges as back up to engineering controls. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

***Protection of hands:**

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Always ensure that gloves are inspected before use.

Selection of protective gloves must include consideration of the penetration times along with rates of diffusion and degradation. The selected glove should comply with the specifications of EU Directive 89/686/EEC and the standard EN374 derived from it.

(Contd. on page 4)

Safety data sheet
according to 1907/2006/EC, Article 31

Printing 27.01.2020

version number 5

Revision: 27.01.2020

Trade name: S3 KINEMATIC VISCOSITY,

(Contd. of page 3)

***Material of gloves**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.

***Penetration time of glove material**

In the absence of data above, the exact break through time has to be sourced from the manufacturer of the protective gloves and has to be observed.

***Eye protection:**

Goggles recommended during refilling: Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU)

***Engineering Controls:** No further data; see item 7.

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties**General Information*****Appearance:**

Form:	Liquid
Colour:	Colourless
*Odour:	Characteristic
*Odour threshold:	Not determined.

***pH-value:** Not determined.

***Change in condition**

Melting point/freezing point:	0 °C
Initial boiling point and boiling range:	Undetermined.

***Flash point:** Not applicable.

***Flammability (solid, gas):** Not applicable.

***Ignition temperature:** 200 °C

***Decomposition temperature:** Not determined.

***Auto-ignition temperature:** Product is not selfigniting.

***Explosive properties:** Product does not present an explosion hazard.

***Explosion limits:**

Lower:	0.6 Vol %
Upper:	40 Vol %

***Vapour pressure at 20 °C:** 23 hPa

***Density:** Not determined.

***Relative density** Not determined.

***Vapour density** Not determined.

***Evaporation rate** Not determined.

***Solubility/ miscibility with water:**

Not miscible or difficult to mix.

***Partition coefficient: n-octanol/water:** Not determined.

***Viscosity:**

Dynamic:	Not determined.
Kinematic:	Not determined.

(Contd. on page 5)

Safety data sheet
according to 1907/2006/EC, Article 31

Printing 27.01.2020

version number 5

Revision: 27.01.2020

Trade name: S3 KINEMATIC VISCOSITY,

(Contd. of page 4)

***9.2 Other information**

No further relevant information available.

SECTION 10: Stability and reactivity

- *10.1 **Reactivity** No further relevant information available.
- *10.2 **Chemical stability**
- ***Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- *10.3 **Possibility of hazardous reactions** No dangerous reactions known.
- *10.4 **Conditions to avoid** No further relevant information available.
- *10.5 **Incompatible materials:** No further relevant information available.
- *10.6 **Hazardous decomposition products:** No dangerous decomposition products known.

SECTION 11: Toxicological information

- *11.1 **Information on toxicological effects**
- ***Acute toxicity** Based on available data, the classification criteria are not met.
- ***Primary irritant effect:**
- ***Skin corrosion/irritation** Based on available data, the classification criteria are not met.
- ***Serious eye damage/irritation** Based on available data, the classification criteria are not met.
- ***Respiratory or skin sensitisation** Based on available data, the classification criteria are not met.
- ***CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)**
- ***Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- ***Carcinogenicity** Based on available data, the classification criteria are not met.
- ***Reproductive toxicity** Based on available data, the classification criteria are not met.
- ***STOT-single exposure** Based on available data, the classification criteria are not met.
- ***STOT-repeated exposure** Based on available data, the classification criteria are not met.
- ***Aspiration hazard**
May be fatal if swallowed and enters airways.

*** SECTION 12: Ecological information**

- *12.1 **Toxicity**
- ***Aquatic toxicity:** No further relevant information available.
- *12.2 **Persistence and degradability** No further relevant information available.
- *12.3 **Bioaccumulative potential** No further relevant information available.
- *12.4 **Mobility in soil** No further relevant information available.
- ***Additional ecological information:**
- ***General notes:**
Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
- *12.5 **Results of PBT and vPvB assessment**
- ***PBT:** Not applicable.
- ***vPvB:** Not applicable.
- *12.6 **Other adverse effects** No further relevant information available.

SECTION 13: Disposal considerations

- *13.1 **Waste treatment methods**
- ***Recommendation**
Must not be disposed together with household garbage. Do not allow product to reach sewage system.

(Contd. on page 6)

IE

Safety data sheet
according to 1907/2006/EC, Article 31

Printing 27.01.2020

version number 5

Revision: 27.01.2020

Trade name: S3 KINEMATIC VISCOSITY,

(Contd. of page 5)

Uncleaned packaging:**Recommendation:** Disposal must be made according to official regulations.**SECTION 14: Transport information*****14.1 UN-Number*****ADR, IMDG, IATA** Void***14.2 UN proper shipping name*****ADR, IMDG, IATA** Void***14.3 Transport hazard class(es)*****ADR, IMDG, IATA*****Class** Void***14.4 Packing group*****ADR, IMDG, IATA** Void***14.5 Environmental hazards:*****Marine pollutant:** No***14.6 Special precautions for user**

Not applicable.

***14.7 Transport in bulk according to Annex II of Marpol and the IBC Code**

Not applicable.

***Transport/Additional information:**

Not dangerous according to the above specifications.

***UN "Model Regulation":**

Void

SECTION 15: Regulatory information***15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture*****Directive 2012/18/EU*****Named dangerous substances - ANNEX I** None of the ingredients is listed.***REGULATION (EC) No 1907/2006 ANNEX XVII** Conditions of restriction: 3***National regulations:*****Waterhazard class:** Water hazard class 1 (Self-assessment): slightly hazardous for water.***15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.**SECTION 16: Other information**

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

***Relevant Phrases:**

H304 May be fatal if swallowed and enters airways.

Department issuing SDS:** Health and SafetyContact:** sds@reagecon.ie***Abbreviations and acronyms:**

REACH (Registration, Evaluation, Authorisation and restriction of Chemicals)

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

(Contd. on page 7)

Safety data sheet
according to 1907/2006/EC, Article 31

Printing 27.01.2020

version number 5

Revision: 27.01.2020

Trade name: S3 KINEMATIC VISCOSITY,

CAS: Chemical Abstracts Service (division of the American Chemical Society)
PBT: Persistent, Bioaccumulative and Toxic
vPvB: very Persistent and very Bioaccumulative
Asp. Tox. 1: Aspiration hazard – Category 1

***NO**

(Contd. of page 6)

IE