

Safety data sheet
according to 1907/2006/EC, Article 31

Printing 07.02.2024

version number 5 (replaces version 4)

Revision: 07.02.2024

SECTION 1: Identification of the substance/mixture and of the company/undertaking***1.1 Product identifier**

***Trade name:** S3 KINEMATIC VISCOSITY,
DYNAMIC VISCOSITY AND DENSITY STANDARD TO 100°C

***Article number:** REVISC-S3

***UFI:** 8HU2-Q00T-F00S-TK75

***1.2 Relevant identified uses of the substance or mixture and uses advised against**

No further relevant information available.

***Application of the substance / the mixture** Laboratory Chemicals

1.3 Details of the supplier of the safety data sheet**Manufacturer/Supplier:**

Reagecon Diagnostics Ltd.

Shannon Free Zone,

Shannon,

Co. Clare,

Ireland.

Tel +353 61 472622

Fax +353 61 472642

***Further information obtainable by contacting:** sds@reagecon.ie

***1.4 Emergency telephone number:**

National Poisons Information Centre: +353 (1) 809 2166 (8.00 a.m. to 10.00 p.m. 7 days a week)

Healthcare Professionals: +353 (1) 809 2566 (24 hour service)

For Hazardous Materials [or Dangerous Goods] Incident

Spill, Leak, Fire, Exposure, or Accident

Call CHEMTREC

For Ireland call +(353)-19014670

For Outside Ireland call +1 703-741-5970 / 1-800-424-9300 CCN849800

SECTION 2: Hazards identification***2.1 Classification of the substance or mixture**

***Classification according to Regulation (EC) No 1272/2008**



GHS08 health hazard

Asp. Tox. 1 H304 May be fatal if swallowed and enters airways.

***2.2 Label elements**

***Labelling according to Regulation (EC) No 1272/2008**

The product is classified and labelled according to the CLP regulation.

***Hazard pictograms**

GHS08

***Signal word** Danger

***Hazard-determining components of labelling:**

Distillates (petroleum), hydrotreated light
dodecane

(Contd. on page 2)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing 07.02.2024

version number 5 (replaces version 4)

Revision: 07.02.2024

**Trade name: S3 KINEMATIC VISCOSITY,
DYNAMIC VISCOSITY AND DENSITY STANDARD TO 100°C**

(Contd. of page 1)

***Hazard statements**

H304 May be fatal if swallowed and enters airways.

***Precautionary statements**

P301+P310 IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

P331 Do NOT induce vomiting.



P405 Store locked up.

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

2.3 Other hazards**Results of PBT and vPvB assessment*****PBT:** Not applicable.***vPvB:** Not applicable.

SECTION 3: Composition/information on ingredients

3.2 Mixtures**Description:** Mixture of substances listed below with nonhazardous additions.***Dangerous components:**

EC number: 927-632-8	Distillates (petroleum), hydrotreated light	 Asp. Tox. 1, H304	50-100%
CAS: 112-40-3 EINECS: 203-967-9	dodecane	 Asp. Tox. 1, H304	10-25%

***Additional information:** For the wording of the listed hazard phrases refer to section 16.

SECTION 4: First aid measures

4.1 Description of first aid measures**After inhalation:**

Provide fresh air, warmth and rest. When breathing is difficult, properly trained personnel may assist affected person by administering oxygen. Obtain medical attention if any discomfort continues.

After skin contact:** Generally the product does not irritate the skin.After eye contact:**

Promptly wash eyes with plenty of water for up to 15 minutes. Open eyes wide apart and rinse well to remove any contact lenses. Do not remove contact lenses by hand. Continue to rinse. Get medical attention if symptoms persist.

***After swallowing:**

Do not induce vomiting; call for medical help immediately. Rinse mouth thoroughly with water and give large amounts of water to drink. Never give anything by mouth to an unconscious person.

Caution if victim vomits. Risk of aspiration. Keep airways free. Call a physician immediately.

4.2 Most important symptoms and effects, both acute and delayed** No further relevant information available.4.3 Indication of any immediate medical attention and special treatment needed**

No further relevant information available.

SECTION 5: Firefighting measures

5.1 Extinguishing media**Suitable extinguishing agents:**

Indications shall be given whether any extinguishing media are inappropriate for a particular situation involving the substance or mixture

Use fire extinguishing methods suitable to surrounding conditions.

(Contd. on page 3)

Safety data sheet
according to 1907/2006/EC, Article 31

Printing 07.02.2024

version number 5 (replaces version 4)

Revision: 07.02.2024

**Trade name: S3 KINEMATIC VISCOSITY,
DYNAMIC VISCOSITY AND DENSITY STANDARD TO 100°C**

(Contd. of page 2)

***5.2 Special hazards arising from the substance or mixture** No further relevant information available.

***5.3 Advice for firefighters**

***Protective equipment:**

Fire fighters should wear complete protective clothing including self-contained breathing apparatus. Do not breathe fumes. Keep containers cool by spraying with water if exposed to fire. Avoid run off to waterways and sewers

SECTION 6: Accidental release measures

***6.1 Personal precautions, protective equipment and emergency procedures**

Wear protective equipment as described in Section 8 below. Keep unprotected persons away.

***6.2 Environmental precautions:** Do not allow to enter sewers/ surface or ground water.

***6.3 Methods and material for containment and cleaning up:**

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

Dispose contaminated material as waste according to section 13.

Ensure adequate ventilation.

***6.4 Reference to other sections**

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

SECTION 7: Handling and storage

***7.1 Precautions for safe handling**

Ensure good ventilation/exhaustion at the workplace.

Prevent formation of aerosols.

***Information about fire - and explosion protection:** No special measures required.

***7.2 Conditions for safe storage, including any incompatibilities**

***Storage:**

***Requirements to be met by storerooms and receptacles:** No special requirements.

***Information about storage in one common storage facility:** Not required.

***Further information about storage conditions:** None.

***7.3 Specific end use(s)** No further relevant information available.

SECTION 8: Exposure controls/personal protection

***8.1 Control parameters**

***Ingredients with limit values that require monitoring at the workplace:**

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

***Additional information:** The lists valid during the making were used as basis.

***8.2 Exposure controls**

***Appropriate engineering controls** No further data; see section 7.

***Individual protection measures, such as personal protective equipment**

***General protective and hygienic measures:**

Keep away from foodstuffs, beverages and feed.

Wash hands before breaks and at the end of work.

(Contd. on page 4)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing 07.02.2024

version number 5 (replaces version 4)

Revision: 07.02.2024

**Trade name: S3 KINEMATIC VISCOSITY,
DYNAMIC VISCOSITY AND DENSITY STANDARD TO 100°C**

(Contd. of page 3)

***Respiratory protection:**

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.

Where risk assessment shows air-purifying respirators are appropriate use a respirator with multi-purpose combination (US) or type ABEK (EN14387) respirator cartridges as back up to engineering controls. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

***Hand protection**

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Always ensure that gloves are inspected before use.

Selection of protective gloves must include consideration of the penetration times along with rates of diffusion and degradation. The selected glove should comply with the specifications of EU Directive 89/686/EEC and the standard EN374 derived from it.

***Material of gloves**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.

***Penetration time of glove material**

In the absence of data above, the exact break through time has to be sourced from the manufacturer of the protective gloves and has to be observed.

***Eye/face protection**

Goggles recommended during refilling: Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU)

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties**General Information*****Physical state**

Liquid

***Colour:**

Colourless

***Odour:**

Characteristic

***Odour threshold:**

Not determined.

***Melting point/freezing point:**

0 °C

***Boiling point or initial boiling point and boiling range**

Undetermined.

***Flammability**

Not applicable.

Lower and upper explosion limit**Lower:**

0.6 Vol %

***Upper:**

40 Vol %

***Flash point:**

Not applicable.

***Auto-ignition temperature:**

200 °C

***Decomposition temperature:**

Not determined.

***pH**

Not determined.

Viscosity:**Kinematic viscosity**

Not determined.

***Dynamic:**

Not determined.

Solubility**water:**

Not miscible or difficult to mix.

***Partition coefficient n-octanol/water (log value)**

Not determined.

***Vapour pressure at 20 °C:**

23 hPa

(Contd. on page 5)

Safety data sheet
according to 1907/2006/EC, Article 31

Printing 07.02.2024

version number 5 (replaces version 4)

Revision: 07.02.2024

**Trade name: S3 KINEMATIC VISCOSITY,
DYNAMIC VISCOSITY AND DENSITY STANDARD TO 100°C**

(Contd. of page 4)

***Density and/or relative density**

*Density:	<i>Not determined.</i>
*Relative density	<i>Not determined.</i>
*Vapour density	<i>Not determined.</i>

***9.2 Other information**

*Appearance:	
*Form:	<i>Liquid</i>
*Important information on protection of health and environment, and on safety.	
*Ignition temperature:	<i>Product is not selfigniting.</i>
*Explosive properties:	<i>Product does not present an explosion hazard.</i>
*Change in condition	
*Evaporation rate	<i>Not determined.</i>

***Information with regard to physical hazard classes**

*Explosives	<i>Void</i>
*Flammable gases	<i>Void</i>
*Aerosols	<i>Void</i>
*Oxidising gases	<i>Void</i>
*Gases under pressure	<i>Void</i>
*Flammable liquids	<i>Void</i>
*Flammable solids	<i>Void</i>
*Self-reactive substances and mixtures	<i>Void</i>
*Pyrophoric liquids	<i>Void</i>
*Pyrophoric solids	<i>Void</i>
*Self-heating substances and mixtures	<i>Void</i>
*Substances and mixtures, which emit flammable gases in contact with water	<i>Void</i>
*Oxidising liquids	<i>Void</i>
*Oxidising solids	<i>Void</i>
*Organic peroxides	<i>Void</i>
*Corrosive to metals	<i>Void</i>
*Desensitised explosives	<i>Void</i>

SECTION 10: Stability and reactivity

- *10.1 Reactivity** *No further relevant information available.*
- *10.2 Chemical stability**
- *Thermal decomposition / conditions to be avoided:** *No decomposition if used according to specifications.*
- *10.3 Possibility of hazardous reactions** *No dangerous reactions known.*
- *10.4 Conditions to avoid** *No further relevant information available.*
- *10.5 Incompatible materials:** *No further relevant information available.*
- *10.6 Hazardous decomposition products:** *No dangerous decomposition products known.*

SECTION 11: Toxicological information

- *11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008**
- *Acute toxicity** *Based on available data, the classification criteria are not met.*
- *Skin corrosion/irritation** *Based on available data, the classification criteria are not met.*

(Contd. on page 6)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing 07.02.2024

version number 5 (replaces version 4)

Revision: 07.02.2024

**Trade name: S3 KINEMATIC VISCOSITY,
DYNAMIC VISCOSITY AND DENSITY STANDARD TO 100°C**

(Contd. of page 5)

- ***Serious eye damage/irritation** Based on available data, the classification criteria are not met.
- ***Respiratory or skin sensitisation** Based on available data, the classification criteria are not met.
- ***Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- ***Carcinogenicity** Based on available data, the classification criteria are not met.
- ***Reproductive toxicity** Based on available data, the classification criteria are not met.
- ***STOT-single exposure** Based on available data, the classification criteria are not met.
- ***STOT-repeated exposure** Based on available data, the classification criteria are not met.
- ***Aspiration hazard** May be fatal if swallowed and enters airways.

***11.2 Information on other hazards**

***Endocrine disrupting properties**

None of the ingredients is listed.

SECTION 12: Ecological information

- ***12.1 Toxicity**
- ***Aquatic toxicity:** No further relevant information available.
- ***12.2 Persistence and degradability** No further relevant information available.
- ***12.3 Bioaccumulative potential** No further relevant information available.
- ***12.4 Mobility in soil** No further relevant information available.
- ***12.5 Results of PBT and vPvB assessment**
- ***PBT:** Not applicable.
- ***vPvB:** Not applicable.
- ***12.6 Endocrine disrupting properties**
The product does not contain substances with endocrine disrupting properties.
- ***12.7 Other adverse effects**
- ***Additional ecological information:**
- ***General notes:**
Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

SECTION 13: Disposal considerations

- ***13.1 Waste treatment methods**
- ***Recommendation**
Must not be disposed together with household garbage. Do not allow product to reach sewage system.
- ***Uncleaned packaging:**
- ***Recommendation:** Disposal must be made according to official regulations.

SECTION 14: Transport information

- | | |
|---------------------------------------|------|
| * 14.1 UN number or ID number | |
| * ADR, IMDG, IATA | Void |
| * 14.2 UN proper shipping name | |
| * ADR, IMDG, IATA | Void |

(Contd. on page 7)

UL

Safety data sheet
according to 1907/2006/EC, Article 31

Printing 07.02.2024

version number 5 (replaces version 4)

Revision: 07.02.2024

**Trade name: S3 KINEMATIC VISCOSITY,
DYNAMIC VISCOSITY AND DENSITY STANDARD TO 100°C**

(Contd. of page 6)

*14.3 Transport hazard class(es)	
*ADR, IMDG, IATA	
*Class	Void
*14.4 Packing group	
*ADR, IMDG, IATA	Void
*14.5 Environmental hazards:	
*Marine pollutant:	No
*14.6 Special precautions for user	Not applicable.
*14.7 Maritime transport in bulk according to IMO instruments	Not applicable.
*Transport/Additional information:	Not dangerous according to the above specifications.
*UN "Model Regulation":	Void

SECTION 15: Regulatory information

*15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

*Directive 2012/18/EU

*Named dangerous substances - ANNEX I None of the ingredients is listed.

*REGULATION (EC) No 1907/2006 ANNEX XVII Conditions of restriction: 3

*DIRECTIVE 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment – Annex II

None of the ingredients is listed.

*REGULATION (EU) 2019/1148

*Annex I - RESTRICTED EXPLOSIVES PRECURSORS (Upper limit value for the purpose of licensing under Article 5(3))

None of the ingredients is listed.

*Annex II - REPORTABLE EXPLOSIVES PRECURSORS

None of the ingredients is listed.

*Regulation (EC) No 273/2004 on drug precursors

None of the ingredients is listed.

*Regulation (EC) No 111/2005 laying down rules for the monitoring of trade between the Community and third countries in drug precursors

None of the ingredients is listed.

*National regulations:

*Waterhazard class: Water hazard class 1 (Self-assessment): slightly hazardous for water.

*15.2 Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

(Contd. on page 8)

UL

Safety data sheet
according to 1907/2006/EC, Article 31

Printing 07.02.2024

version number 5 (replaces version 4)

Revision: 07.02.2024

**Trade name: S3 KINEMATIC VISCOSITY,
DYNAMIC VISCOSITY AND DENSITY STANDARD TO 100°C**

(Contd. of page 7)

***Relevant Phrases:**

H304 May be fatal if swallowed and enters airways.

Department issuing SDS:** Health and SafetyContact:** *sds@reagecon.ie****Version number of previous version:** 4***Abbreviations and acronyms:**

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

Asp. Tox. 1: Aspiration hazard – Category 1

UL