

Safety data sheet

according to Regulation (EC) No 1907/2006 (REACH)

Printing date 20.06.2017

Version number 3

Revision: 19.06.2017

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

- Trade name: **ISA/Ag/FI**
- Article number: 140140

Description: ISA solution for the determination of silver in in fixer solution

1.2 Relevant identified uses of the substance or mixture and uses advised against:

- Product category: PC21 Laboratory chemicals
- Process category: PROC15 Use as laboratory reagent
- Application of the substance / the preparation: Sample conditioning solution

1.3 Details of the supplier of the safety data sheet

Manufacturer/Supplier:

Xylem Analytics Germany GmbH
WTW
Dr.-Karl-Slevogt-Straße 1
82362 Weilheim
Germany
Tel. +49 881 183-0

Further information obtainable from: E-mail: Info.WTW@Xyleminc.com

1.4 Emergency telephone number: Chemtrec: (USA & Canada) 800-424-9300 (International) 001 703-527-3887

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

Classification according to Regulation (EC) No 1272/2008: The product is not classified according to the CLP regulation.

2.2 Label elements:

- Labelling according to Regulation (EC) No 1272/2008: Void
 - Hazard pictograms: Void
 - Signal word: Void
 - Hazard statements: Void

2.3 Other hazards No further relevant information available.

SECTION 3: Composition/information on ingredients

3.2 Mixture

Description:

Mixture of substances listed below with nonhazardous additions.
Water, sodium sulfite

Dangerous components: Void

Additional information: For the wording of the listed hazard phrases refer to section 16.

SECTION 4: First aid measures

4.1 Description of first aid measures

- After inhalation: Supply fresh air.
- After skin contact: Rinse skin with water/shower.
- After eye contact: Rinse opened eye for several minutes under running water.
- After swallowing:
Rinse out mouth and then drink plenty of water.
Get medical advice/attention if you feel unwell.

4.2 Most important symptoms and effects, both acute and delayed No further relevant information available.

(Contd. on page 2)

Safety data sheet

according to Regulation (EC) No 1907/2006 (REACH)

Printing date 20.06.2017

Version number 3

Revision: 19.06.2017

Trade name: ISA/Ag/FI

(Contd. of page 1)

- **4.3 Indication of any immediate medical attention and special treatment needed** No further relevant information available.

SECTION 5: Firefighting measures

- **5.1 Extinguishing media**
 - **Suitable extinguishing agents:** The product is not flammable. Extinguishing agent to suit environment.
- **5.2 Special hazards arising from the substance or mixture**
In case of fire, the following can be released:
Sulphur oxides (SOx)
- **5.3 Advice for firefighters**
 - **Protective equipment:** Wear self-contained respiratory protective device.

SECTION 6: Accidental release measures

- **6.1 Personal precautions, protective equipment and emergency procedures**
No special measurements required when handled properly.
- **6.2 Environmental precautions:**
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
- **6.3 Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
Wash off residuals with water.
- **6.4 Reference to other sections** No dangerous substances are released.

SECTION 7: Handling and storage

- **7.1 Precautions for safe handling**
No special precautions are necessary if used correctly.
Do not leave the container open when not used.
 - **Information about fire - and explosion protection:** No special measures required.
- **7.2 Conditions for safe storage, including any incompatibilities**
 - **Requirements to be met by storerooms and receptacles:** No special requirements.
 - **Information about storage in one common storage facility:** Not required.
 - **Further information about storage conditions:** Store tightly sealed at temperatures between 15 °C and 25 °C.
- **7.3 Specific end use(s)** No further relevant information available.

SECTION 8: Exposure controls/personal protection

- **Additional information about design of technical facilities:** No further data; see item 7.
- **8.1 Control parameters**
 - **Ingredients with limit values that require monitoring at the workplace:**
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
 - **Additional information:** The lists valid during the making were used as basis.
- **8.2 Exposure controls**
 - **Personal protective equipment:**
 - **General protective and hygienic measures:**
Avoid contact with the eyes and skin.
Do not inhale gases / fumes / aerosols.
Take off contaminated clothing and wash it before reuse.
Wash hands before breaks and at the end of work.
 - **Respiratory protection:** Not required.
 - **Protection of hands:** Protective gloves recommended for frequent or long-term skin contact.
 - **Material of gloves** Nitrile rubber, NBR

(Contd. on page 3)

Safety data sheet

according to Regulation (EC) No 1907/2006 (REACH)

Printing date 20.06.2017

Version number 3

Revision: 19.06.2017

Trade name: ISA/Ag/FI

(Contd. of page 2)

· Eye protection: Safety glasses

SECTION 9: Physical and chemical properties

· 9.1 Information on basic physical and chemical properties	
· General Information	
· Appearance:	
· Form:	Fluid
· Colour:	Colourless
· Odour:	Odourless
· pH-value at 20 °C:	6.25
· Change in condition	
· Melting point/freezing point:	Undetermined.
· Initial boiling point and boiling range:	100 °C
· Flash point:	Not applicable.
· Explosive properties:	Product does not present an explosion hazard.
· Vapour pressure at 20 °C:	23 hPa
· Density at 20 °C:	1.07 g/cm ³
· Solubility in / Miscibility with	
· water:	Fully miscible.
· Viscosity:	
· Dynamic:	Not determined.
· Kinematic:	Not determined.
· 9.2 Other information	No further relevant information available.

SECTION 10: Stability and reactivity

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability**
 - **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **10.3 Possibility of hazardous reactions** No dangerous reactions known.
- **10.4 Conditions to avoid** No further relevant information available.
- **10.5 Incompatible materials:**
 - Acids
 - Oxidizing agents
 - Atmospheric oxygen
- **10.6 Hazardous decomposition products:** In case of fire, see section 5.

SECTION 11: Toxicological information

- **11.1 Information on toxicological effects**
 - **Acute toxicity** Based on available data, the classification criteria are not met.
 - **Primary irritant effect:**
 - **Skin corrosion/irritation** Local irritation possible.
 - **Serious eye damage/irritation** Irritation possible.
 - **Respiratory or skin sensitisation** Based on available data, the classification criteria are not met.
 - **CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)**
 - **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
 - **Carcinogenicity** Based on available data, the classification criteria are not met.
 - **Reproductive toxicity** Based on available data, the classification criteria are not met.
 - **STOT-single exposure** Based on available data, the classification criteria are not met.
 - **STOT-repeated exposure** Based on available data, the classification criteria are not met.

(Contd. on page 4)

Safety data sheet

according to Regulation (EC) No 1907/2006 (REACH)

Printing date 20.06.2017

Version number 3

Revision: 19.06.2017

Trade name: ISA/Ag/FI

(Contd. of page 3)

- **Aspiration hazard** Based on available data, the classification criteria are not met.

SECTION 12: Ecological information

- **12.1 Toxicity**
 - **Aquatic toxicity:** No further relevant information available.
- **12.2 Persistence and degradability** No further relevant information available.
- **12.3 Bioaccumulative potential** No further relevant information available.
- **12.4 Mobility in soil** No further relevant information available.
 - **Ecotoxicological effects:**
 - **Remark:** No environmental hazards are to be expected when the product is handled properly.
 - **Additional ecological information:**
 - **General notes:**
 - Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water
 - Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
 - **Additional ecological information:**
 - Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water
 - Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
- **12.5 Results of PBT and vPvB assessment** Not applicable.
- **12.6 Other adverse effects**
 - **General notes:**
 - Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water
 - Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

SECTION 13: Disposal considerations

- **13.1 Waste treatment methods**
 - **Recommendation**
 - Disposal must comply with the relevant local regulations. We recommend that you contact either the authorities in charge or approved waste disposal companies which will advise you on how to dispose the special waste.
 - **Uncleaned packaging:**
 - **Recommendation:**
 - Empty contaminated packagings thoroughly. They may be recycled after thorough and proper cleaning.
 - **Recommended cleansing agents:** Water, if necessary together with cleansing agents.

SECTION 14: Transport information

· 14.1 UN-Number	
· ADR/RID, ADN, IMDG, IATA	Void
· 14.2 UN proper shipping name	
· ADR/RID, ADN, IMDG, IATA	Void
· 14.3 Transport hazard class(es)	
· ADR/RID, ADN, IMDG, IATA	
· Class	Void
· 14.4 Packing group	
· ADR/RID, IMDG, IATA	Void
· 14.5 Environmental hazards:	Not applicable.
· 14.6 Special precautions for user	Not applicable.
· 14.7 Transport in bulk according to Annex II of Marpol and the IBC Code	Not applicable.
· UN "Model Regulation":	Void

EU
(Contd. on page 5)

Safety data sheet

according to Regulation (EC) No 1907/2006 (REACH)

Printing date 20.06.2017

Version number 3

Revision: 19.06.2017

Trade name: ISA/Ag/FI

(Contd. of page 4)

* SECTION 15: Regulatory information

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
 - **Directive 2012/18/EU**
 - **Named dangerous substances - ANNEX I** None of the ingredients is listed.
- **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

SECTION 16: Other information

- **Abbreviations and acronyms:**

RID: Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail)

IATA-DGR: Dangerous Goods Regulations by the "International Air Transport Association" (IATA)

ICAO: International Civil Aviation Organisation

ICAO-TI: Technical Instructions by the "International Civil Aviation Organisation" (ICAO)

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

- *** Data compared to the previous version altered.**